



## ProSwift® SAX-1S, WAX-1S, SCX-1S, and WCX-1S QUICKSTART

ProSwift Monolith columns offer superior separation of large biomolecules. Conditioning of the column bed is **required** prior to initial use and after long-term storage. The QuickStart process will ensure extended column lifetime and reproducibility.

### I. Preparation:

#### a. Eluent Preparation

The following eluents are recommended, but the column may be used with any eluents suitable for analysis. Typically Eluent A is a low salt buffer and Eluent B is a high salt buffer.

<u>Eluent</u>	<u>Anion Exchange</u>	<u>Cation Exchange</u>
A	0.01 M Tris·HCl buffer pH 7.6	0.01M Na <sub>2</sub> PO <sub>4</sub> buffer, pH 7.6
B	1 M NaCl in 0.01 M Tris •HCl buffer pH 7.6	1 M NaCl in 0.01 M Na <sub>2</sub> PO <sub>4</sub> buffer, pH 7.6

#### b. Column Installation

Install the column on the instrument in the correct flow direction.



***Sudden increases in flow rates may damage monolithic columns. Always increase the flow rate slowly using a linear flow gradient or stepwise increments in flow rate.***

***If the eluent composition generates back pressure in excess of the maximum pressure, reduce the flow rate to ensure the upstream back pressure is less than the maximum pressure.***

***The maximum pressure limit for 4.6mm is 1000 psi (6.9 MPa) and for 1mm is 2000 psi (13.8 MPa).***

### II. Flow Rate Start-Up (Ramping)

Using a linear or stepwise flow gradient, increase the flow rate of Eluent A starting from 0.00 mL/min to the desired flow rate using the flow rates given below.

- 4.6mm Column: Use a rate of  $\leq 0.50$  mL/min, every 30 seconds
- 1mm Column: Use a rate of  $\leq 0.10$  mL/min every minute.

### III. Column Conditioning - Use the guidelines below to determine the proper startup conditions:

#### 4.6 mm:

##### Removal of Storage Solution:

- Using the desired flow rate, run a 8 mL (~10 column volume) binary gradient from 100% A to 100% B.
- Pump another 8 mL (~10 column volumes) of 100% B through the column.

##### Column Equilibration:

- Equilibration from 100% B to the starting eluent composition should include at least a 1 minute reverse gradient to the initial conditions.
- Pump at least 8 mL (~10 column volumes) of this eluent composition through the column.

#### 1 mm:

##### Removal of Storage Solution:

- Using the desired flow rate, run a 0.4 mL (~10 column volume) binary gradient from 100% A to 100% B.
- Pump another 0.4 mL (~10 column volumes) 100% B through the column.

##### Column Equilibration:

- Equilibration from 100% B to the starting eluent composition should include at least a 1 minute reverse gradient to the initial conditions.
- Pump at least 0.4 mL (~10 column volumes) of this eluent composition through the column.

### IV. Storage:

- For short-term storage, <3 days, store the column in the initial buffer compositions.
- For long-term storage, >3 days, use the following solution to avoid microbial growth on the column.

##### Anion Exchange

0.05 M NaCl in 0.01 M Tris·HCl, pH 7.6 + 0.1% NaN<sub>3</sub>

##### Cation Exchange

0.05 M NaCl in 0.01 M Na<sub>2</sub>PO<sub>4</sub>, pH 7.6 + 0.1% NaN<sub>3</sub>

- In all cases the column should be tightly sealed with end plugs to prevent the column from drying out.

***For additional information, please refer to the manual, ProSwift Ion Exchange Product Manual Doc. No. 065122; see Section "System Requirements" regarding the use of stainless steel HPLC systems.***

Now sold under the  
Thermo Scientific brand

**Thermo**  
SCIENTIFIC