

METHOD 3810

HEADSPACE

1.0 SCOPE AND APPLICATION

1.1 Method 3810 was formerly Method 5020 in the second edition of this manual.

1.2 Method 3810 is a static headspace technique for extracting volatile organic compounds from samples. It is a simple method that allows large numbers of samples to be screened in a relatively short period of time. It is ideal for screening samples prior to using the purge-and-trap method. Detection limits for this method may vary widely among samples because of the large variability and complicated matrices of waste samples. The method works best for compounds with boiling points of less than 125°C. The sensitivity of this method will depend on the equilibria of the various compounds between the vapor and dissolved phases.

1.3 Due to the variability of this method, this procedure is recommended for use only as a screening procedure for other, more accurate determinative methods (Methods 8010, 8015, 8020, 8030, and 8240).

2.0 SUMMARY OF METHOD

2.1 The sample is collected in sealed glass containers and allowed to equilibrate at 90°C. A sample of the headspace gas is withdrawn with a gas-tight syringe for screening analysis using the conditions specified in one of the GC or GC/MS determinative methods (8010, 8015, 8020, 8030, or 8240).

3.0 INTERFERENCES

3.1 Samples can be contaminated by diffusion of volatile organics (particularly chlorofluorocarbons and methylene chloride) through the sample container septum during shipment and storage. A field sample blank prepared from reagent water and carried through sampling and subsequent storage and handling can serve as a check on such contamination.

3.2 Contamination by carryover can occur whenever high-level and low-level samples are sequentially analyzed. To reduce carryover, the sample syringe must be rinsed out between samples with reagent water. Whenever an unusually concentrated sample is encountered, it should be followed by an analysis of reagent water. It may be necessary to wash out the syringe with detergent, rinse with distilled water, and dry in a 105°C oven between analyses.

3.3 Before processing any samples, the analyst should demonstrate daily through the analysis of an organic-free water or solvent blank that the entire analytical system is interference-free.

4.0 APPARATUS AND MATERIALS

4.1 Refer to the specific determinative method for appropriate apparatus and materials.

4.2 Vials: 125-mL Hypo-Vials (Pierce Chemical Co., #12995, or equivalent), four each.

4.3 Septa: Tuf-Bond (Pierce #12720 or equivalent).

4.4 Seals: Aluminum (Pierce #132141 or equivalent).

4.5 Crimper: Hand (Pierce #13212 or equivalent).

4.6 Syringe: 5-mL, gas-tight with shutoff valve and chromatographic needles.

4.7 Microsyringe: 250- or 500-uL.

4.8 Water bath: Heated, with concentric ring cover, capable of temperature control ($\pm 5^{\circ}\text{C}$). The bath should be used in a hood.

5.0 REAGENTS

5.1 Refer to the specific determinative method and Method 8000 for preparation of calibration standards.

6.0 SAMPLE COLLECTION, PRESERVATION, AND HANDLING

6.1 Refer to the introductory material to this chapter, Organic Analytes, Section 4.1.

7.0 PROCEDURE

7.1 Gas chromatographic conditions and Calibration: Refer to the specific determinative method for GC operating conditions and to Method 8000, Section 7.4, for calibration procedures.

7.2 Sample preparation:

7.2.1 Place 10.0 g of a well-mixed waste sample into each of two separate 125-mL septum-seal vials.

7.2.2 Dose one sample vial through the septum with 200 uL of a 50 ng/uL calibration standard containing the compounds of interest. Label this "1-ppm spike."

7.2.3 Dose a separate (empty) 125-mL septum seal vial with 200 uL of the same 50 ng/uL calibration standard. Label this "1-ppm standard."

7.2.4 Place the sample, 1-ppm-spike, and 1-ppm-standard vials into a 90°C water bath for 1 hr. Store the remaining sample vial at 4.0°C for possible future analysis.

7.3 Sample analysis:

7.3.1 While maintaining the vials at 90°C, withdraw 2 mL of the headspace gas with a gas-tight syringe and analyze by direct injection into a GC. The GC should be operated using the same GC conditions listed in the method being screened (8010, 8015, 8020, 8030, or 8240).

7.3.2 Analyze the 1-ppm standard and adjust instrument sensitivity to give a minimum response of at least 2 times the background. Record retention times (RT) and peak areas of compounds of interest.

7.3.3 Analyze the 1-ppm spiked sample in the same manner. Record RTs and peak areas.

7.3.4 Analyze the undosed sample as in Paragraph 7.3.3.

7.3.5 Use the results obtained to determine if the sample requires dilution or methanolic extraction as indicated in Method 5030.

8.0 QUALITY CONTROL

8.1 Before processing any samples, the analyst should demonstrate through the analysis of a distilled water method blank that all glassware and reagents are interference-free. Each time a set of samples is extracted or there is a change in reagents, a method blank should be processed as a safeguard against chronic laboratory contamination. The blank samples should be carried through all stages of the sample preparation and measurement.

8.2 Standard quality assurance practices should be used with this method. Fortified samples should be carried through all stages of sample preparation and measurement; they should be analyzed to validate the sensitivity and accuracy of the analysis. If the fortified waste samples do not indicate sufficient sensitivity to detect less than or equal to 1 ug/g of sample, then the sensitivity of the instrument should be increased.

9.0 METHOD PERFORMANCE

9.1 No data provided.

10.0 REFERENCES

1. Hachenberg, H. and A. Schmidt, Gas Chromatographic Headspace Analysis, Philadelphia: Hayden & Sons Inc., 1979.
2. Friant, S.L. and I.H. Suffet, "Interactive Effects of Temperature, Salt Concentration and pH on Headspace Analysis for Isolating Volatile Trace Organics in Aqueous Environmental Samples," Anal. Chem. 51, 2167-2172, 1979.

METHOD 3810
HEADSPACE METHOD

